


Expert Pain Killer **e-TRON** IN-1300



Random wave therapy

By anesthetic effect, e-tron is very effective on chronic pain relief

After image effect

Strong persistence continues to exercise and treat after the therapy


Anesthetic effect

Anesthetize applied area to reduce the pain and stress from associated treatment such as injection therapy

Various output modes

Three different output mode available (massage, continuous, variety)

e-tron, IN-1300, is highly effective device for people who have chronic muscle fatigue, lum-bago, musculoskeletal system disorder, ar-thralgia, myalgia, discomfort shoulder includ-ing the frozen shoulder



Features of goods

- The function of the afterimage remains: The amazing waveform composition is very effective for chronic pain due to anesthetic effect.
- Last setting save function
- Zero start (safety device) function
- Overcurrent protection device
- 1st class equipment and B-type

Product Constitution

Body 1ea, Power cord 1ea, Suction cup 8ea, Round resin pulp 8ea, Pad line 4ea, Square pad net 8ea, Band (large) 2ea, Band (Small) 4ea, Square resin pulp 8ea, Suction pad line 2ea., Square pad 8ea

Specification

- Rated Voltage and Input Frequency : AC 220V, 60Hz
- Power consumption : 100VA
- Maximum output power : 35mA
- Frequency range : e-Tron Frequency range : Low = 1000~1255Hz, High = High= 1000~1480Hz
- TENS Frequency range : 1~150Hz, 500Hz, 1000Hz
- Vacuum suction cycle : -Vacuum release : 1.3 sec
- Vacuum time : 1.5sec, 2.0sec, 2.5sec, 3.0sec, 4.0sec, 5.0sec
- Auto vacuum time, Continuous vacuum
- Vacuum suction cycle - Auto time : Vacuum release = 1.3sec
- Vacuum time = 5.0sed, 1.0sec
- Mode selectable - 3 Mode : Mod(Intermit), Cont(Continue), Cross
- Usage time : 1 ~ 60 minute
- Dimension / Weight : (W) 516mm x (L) 355mm x (H) 145mm / 9.5Kg

Usage

- Hospital, Clinic, Oriental medicine Clinic, Physiotherapy Clinic..., etc