



ISO 13485:2016



INNOBIZ
Innovation Business Company

MiracleQ

Electro-Magnetic therapy



MiracleQ

Electro-Magnetic therapy

Customized system with best performance and convenience with new pulse magnetic field energy



① Magnetic field applicator

Ergonomic head design provides high energy penetration and intensively stimulates deep pain sources



② Advanced settings for custom treatment

Dial button selection incorporates parameter adjustment settings, Easy Use!



③ Safe environment

Temperature sensor function safety switch, over-current protection, zero-start safety device applied




④ Memory function

Up to 10 memory settings can be stored.



PEMF




PEMF (Pulsed Electro Magnetic Frequency)

Transmission of 'soft but powerful' magnetic pulsed energy into the human body to activate cells deep into muscle and nerve tissue and relieve pain

- Stimulates deep pain sources due to high energy penetration rate
- Implementing a pulse magnetic field of a safe and clean waveform.

1 Pain Relief

Miracle-Q instantly permeates powerful pulsed magnetic fields into the human body to activate cells deep in muscles and nerve tissue and to help ease musculoskeletal disease pain.



2 Body type management

Magnetic field energy stimulates core muscles to help strengthen muscle and maintain body shape. It can be applied to other muscles, and it helps balance the body healthily.



3 MMS (Muscle Magnetic Stimulator)

Miracle-Q helps strengthen the elasticity and repair of pelvic muscle tissue



**Pulse
Magnetic
Stimulation**

**Powerful
7.7 TESLA**

Safety


**Deep
Stimulation**

**Non-
Invasive**


ADVANTAGE

The magnetic field can easily stimulate muscle and nervous tissue without fat and bone interference to the extent that it does not harm the human body.

- Low burden on human body due to non-invasive stimulation
- Short treatment time
- Reduce product efficiency consumables costs
- Applicable to various areas
- Deeper stimulation than conventional electrical stimulation
- Energy transfer without changing clothes



Front - runner of physiotherapy device



Product name (model name)	Miracle Q (IN-3200Z)
Rated voltage / frequency	AC 220V / 30Hz
Power consumption	1500 VA
Output power	7.7 Tesla $\pm$ 20%
Frequency range	1-50Hz
Stimulus pulse width	300 μ s $\pm$ 20%
Stimulus duration	1 ~ 30 sec
Resting between stimuli	1 ~ 20 sec
Dimension(mm)	387(W) X 473(D) X 954(H)/30Kg



394 Nambusunhwan-ro, Yangcheon-gu, Seoul, Republic of Korea (198-5, Sinwol-dong)
T. +82-2-3663-8329 F. +82-2-2699-8328 E. trade@young-inBT.com